

# Interferometric Autocorrelator

## Features

- ◆ Ultrafast Autocorrelation Measurements for 650 - 1100 nm
- ◆ For Pulses from 40 fs to 1 ps or Pulses Down to 15 fs with Precompensation
- ◆ Integrated Alignment Target and Iris for Easy Setup
- ◆ Selectable Gain up to 70 dB to Accommodate Large Range of Input Powers

The FSAC Interferometric Autocorrelator provides approximate measurements of ultrafast pulse widths. It is based on a modified Michelson interferometer with a nonlinear detector at its output.

This autocorrelator is designed for pulse durations ranging from 15 fs to 1 ps (with precompensation), and the scan range can be adjusted from  $\pm 25$  fs to  $\pm 5$  ps. A BNC output connector provides the autocorrelation trace, which can be viewed on any oscilloscope with  $>1.5$  MHz bandwidth. A photodiode amplifier with selectable gain up to 70 dB allows the user to adjust the unit for a large range of input powers.

Intuitive, easy-to-reach controls on the outside of the enclosure, shown in the photo to the right, allow the user to optimize the interferogram. A MATLAB® or Python® script is available for converting interferometric autocorrelation data to an intensity correlation trace; this is ideal for cases where bandwidth is insufficient or the pulse is excessively chirped.



FSAC  
Front View




FSAC  
Back View

## Specifications

Item #	FSAC
Input Wavelength Range	650 - 1100 nm
Input Pulse Duration (FWHM)	40 fs - 1 ps (Without Precompensation) 15 fs - 1 ps (With Precompensation)
Full Scan Range	50 fs - 10 ps ( $\pm 25$ fs to $\pm 5$ ps)
Input Repetition Rate	$>300$ kHz
Noise-Equivalent Sensitivity <sup>a</sup>	$0.1 \text{ W}^2$ at 800 nm for $\varnothing 1$ mm Beam ( $1/e^2$ )
Input Polarization	Horizontal
Input Beam Diameter	$<\varnothing 4$ mm ( $1/e^2$ )
Scan Rate	5 Hz
Internal Dispersion (Nominal)	230 fs <sup>2</sup> at 800 nm (GDD) 345 fs <sup>3</sup> at 800 nm (TOD)

a. Peak Power x Average Power of Input Laser



This interferometric autocorrelation was obtained using the OCTAVIUS-85M-HP. Precompensating mirrors were used to chirp the beam prior to the autocorrelator.