

## Handy Type Model : H350

Angle measurement

Wavefront sensor

6-axis sensor O-PIAS

Laser autocollimator

Laser autocollimator set

Accessories for alignment

Wavelength 408 / 650	WD(mm) ~300	Beam diameter (mm) Φ 1~6	Resolution 8sec.*	Measurement range	±0.5	±1.0	±1.75
*It might be different based on model and beam diameter.						*H350R-C050 + When connected processing unit.	



## ■ Features

Reasonable and compact size.

## ■ High performance and compact size

High quality is condensed in a compact size.  
Working distance is as same as standard model (H499 series).


## ■ Switchable outgoing beam radius

Can be switched the beam radius by pin hole in a compact body.  
 $\varphi 1,2,6$  as standard.  $\varphi 3$  and  $4$  is option.  
※H350R-C175 is for  $\varphi 1,2$  and  $5$


## ■ Equipped with red LD

H350R series equipped with red LD, H350B has blue LD for pick up.

## ■ Application



built-in system (Glass parallelism measurement)



Mirror angle measurement for digital single-lens reflex

High speed·  
High resolution

Technical guide

Selection guide

Compact

Standard

High resolution wide view

Screen

Ghost cancel

2 wavelength

Processing unit

Accessories

	SPEC						
	Red LD-mounted model			Blue LD-mounted model			
Model	H350R-C050	H350R-C100	H350R-C175	H350B-C050	H350B-C100		
Measurement range	±0.5 deg.	±1.0 deg.	±1.75 deg.	±0.5 deg.	±1.0 deg.		
Resolution*	8 sec.	15 sec.	27 sec.	8 sec.	15 sec.		
Measurement wavelength	650±10nm			408±10nm			
Inbuilt light source wavelength	Semiconductor laser 650±10nm (with ON/OFF switch)			Semiconductor laser 408±10nm (with ON/OFF switch)			
Laser output	1.0mW or less (Class 2)						
Outgoing beam diameter	$\varphi 1, 2, 6\text{mm}$ * * : Selected depending on the pin-hole ( $\varphi 3$ , pin-hole for 4mm are sold separately)						
Using The measurement object distance $\varphi 1, 2, 3\text{mm}$	0~300mm	0~200mm	0~120mm	0~300mm	0~200mm		
Using The measurement object distance $\varphi 4, 5, 6\text{mm}$ **	0~200mm	0~160mm		0~200mm	0~160mm		
Output	NTSC method video output, dedicated cable(1.5m) BNC connection						
Weight	0.2kg (Without AC adaptor)						
Electrical power (Included AC adaptor)	Input:AC100-240V 50/60Hz Output:DC6V 1,000mA	Input:AC100-240V 50/60Hz Output:DC12V 500mA	Input:AC100-240V 50/60Hz Output:DC6V 850mA				
Head size (Not include protruding part)	W40×D50×H80mm		W40×D50×H99mm		W40×D50×H85mm		
Standard features	LD switch, Controled LD light level, CCD sensitivity calibration, switching beam diameters, available mounted at 4 any position						
Accessary	Video cable, AC adaptor, Pin-hole (for $\varphi 1, 2\text{mm}$ ), Instruction manua						
Compatible tilt stage	HB10						

\* Pixel resolution at connected to the processing unit. \*\* H350R-C175 is  $\varphi 1, 2, 5, 6\text{mm}$ .

Wavefront sensor

6-axis sensor O-PIAS

Laser autocollimator

Laser autocollimator set

Accessories for alignment

High speed·  
High resolution

Technical guide

Selection guide

Compact

Standard

High resolution wide view

Screen

Ghost cancel


2 wavelength

Processing unit


Accessories

■ External dimensions


■ H350\*-C050, C100



■ H350R-C175



■ Switch · Connector



■ Model

# H350R-C175/F12E

1 2 3 4

1 Inbuilt light source wavelength

Model	Wavelength $\lambda$
R	650nm
B	408nm

2 Measurement range / Resolution<sup>※1</sup>

Model	Standard(Measuring reflection)	Type (External incidence)
C 0 5 0	$\pm 0.50$ deg. 8 sec.	$\pm 1.00$ deg. 15 sec.
C 1 0 0	$\pm 1.00$ deg. 15 sec.	$\pm 2.00$ deg. 30 sec.
C 1 7 5 <sup>※2</sup>	$\pm 1.75$ deg. 27 sec.	$\pm 3.50$ deg. 53 sec.

4 Incident light

Model	Factory setting
Unshown	Standard(Measuring reflection)
E	External incidence measurement

※1 The resolution in the case of the set including the processing unit.

※2 Only H350R

3 Set format

Model	Product name	Laser autocollimator	Processing unit	LCD monitor	Monitor bracket	Tilt stage	Stand	Accesarry
/ F 1 2	Full-set1200	H350*	HIP-1200	VCM-562HIP	HA10-MB	HB10	HA10	H350-AP
/ F 5 5	Full-set550		HIP-550	VCM-562HIP	HA10-MB	HB10	HA10	H350-AP
/ M 1 2	Monitor and HIP set 1200		HIP-1200	VCM-562HIP	—	—	—	—
/ M 5 5	Monitor and HIP set 550		HIP-550	VCM-562HIP	—	—	—	—
/ H 1 2	HIP set 1200		HIP-1200	—	—	—	—	—
/ H 5 5	HIP set 550		HIP-550	—	—	—	—	—
/ D - H	Half set(Display deg.)		—	VCM-562P	HA10-MB	HB10	HA10	H350-AP Degree displayed film
/ M - H	Half set(Display min.)		—	VCM-562P	HA10-MB	HB10	HA10	H350-AP Minute displayed film
/ D	Monitor set(Display deg.)		—	VCM-562P	—	—	—	Degree displayed film
/ M	Monitor set(Display min.)		—	VCM-562P	—	—	—	Minute displayed film
Unshown	Laser autocollimator		—	—	—	—	—	—

# Optical Sensor

Measurable angle and tilt without contact

## Laser Autocollimator

Angle measurement



### Features

■ An angle is measured instantaneously with non-contact and one point measurement.

■ Long working distance. (Max. 350mm)\*


■ High resolution. (0.0008 degree)\*

\*That differs depending on the product.

### <Measurement comparison>


#### Conventional system

##### ① Displacement sensor system





#### Laser autocollimator system

##### Displacement sensor



##### ② Optical autocollimator system



Every tilt of objects is measured for the easy way

(Laser specular reflection is required when measuring)



Laser autocollimator

High speed·  
High resolution

Technical  
guide

Selection  
guide

Compact

Standard

High resolution  
wide view

Screen

Ghost  
cancel

2 wavelength

Processing  
unit

Accessories

## Configuration of Laser Autocollimator

Measurable tilt with only an autocollimator.

### ■ Appearance



※Full set model

### ■ Standard



Laser autocollimator (LAC)  
• Lineup:25models



Processing unit (HIP)  
• Display tilt digitally



LCD  
• 5.6 inch LCD  
• Cable length:1.5m  
• Connect only NTSC



Tilt stage  
• 2-axis tilt mechanism (with lock)



Monitor bracket  
• Pole and monitor are fixed.



Stand  
• Set up any optical sensor.

Angle measurement

Wavefront sensor

6-axis sensor O-PIAS

Laser autocollimator

Laser autocollimator set

Accessories for alignment

High speed·  
High resolution

Technical guide

Selection guide

Compact

Standard

High resolution wide view

Screen

Ghost cancel

2 wavelength

Processing unit

Accessories

## Set models

### ■ Digital measurement



Full set



High rigidity set




HIPset




Monitor&HIPset

### ■ Visual measurement



Half set




Monitor set

With stand (This is cross-compatible across development, preproduction and manufacturing.)

Without stand (Especially corresponds to manufacturing)


# How to Select

## [Before use]



## [Selection guide]

### ① Measurement range



### ② Pixel resolution

Measurement range ( deg.)		± 0.16	± 0.25	± 0.3	± 0.35	± 0.5	± 1.0	± 1.5	± 1.75	± 2.0
Pixel resolution	Processing unit ( sec.)	3	4	5	6	8	15	23	27	30

※ Resolution of half and monitor set is one tenth that of the measurement range.

※ We would like to recommend H900 when you need pixel resolution in less than 3 seconds.

© P.5-017

### ③ Wavelength


Wavelength	Measurement objects
● 650 nm (Red)	Mirror、 Pure glasses
● 405 nm (Blue)	Not available red LD due to refraction coating and IR cutting filter.

### ④ The measurement object distance

Confirmed the measurement object distance due to model type.

#### <Example>



Auto collimator  
Measurement range :1.0 deg.



## Processing Unit : HIP Series


Tilt of work is displayed in numeric value by connecting with laser autocollimator.

High performance reflected light gravity analysis function is improved the measurement accuracy.



**Application**

This processing unit can be supported productivity improvement.



**HIP Comparison**

HIP Type	HIP-550	HIP-1200
Single measurement	○	○
Multi measurement	×	○
Communication speed	30Hz	5Hz
Runout, axis analysis	×	○
Interface	Serial Parallel	○ ○


**SPEC**

	Typical type	High performance type
		HIP-1200
Angular measurement mode	Single measurement	○
	Multi measurement	○
Measurement application	Off-set tilt measurement	○
	3 beam measurement	○
Centroid analysis	Runout measurement	○
	Center of area	○
External in and output	Luminance centroid	○
	Luminance peak	○
	Averaging	○
	I/O for control	○
	Communication interface	○
		○


※HIP-500/HIP-1000 has stopped being produced. Please contact us for maintenance.

■ Dimensional outline drawings



■ HIP-550



■ HIP-1200



■ Switch · Connector · Connection



\* Blue area is included a unit model.

■ Model

# HIP-1200/016-PE

1 2 3 4

1 Processing unit

Model	Type
HIP-550	Typical type
HIP-1200	High performance type

HIP-550



HIP-1200



\*HIP-500/HIP-1000 has stopped being produced. Please contact us for maintenance.

2 Laser autocollimator Measurement range

Model	Factory setting
016	$\pm 0.16^\circ$
020	$\pm 0.20^\circ$
025	$\pm 0.25^\circ$
030	$\pm 0.30^\circ$
035	$\pm 0.35^\circ$
050	$\pm 0.50^\circ$
100	$\pm 1.00^\circ$
150	$\pm 1.50^\circ$
175	$\pm 1.75^\circ$
200	$\pm 2.00^\circ$

3 Laser autocollimator type

Model	Type
Unshown	H400, H600
P	H350
V	H450

4 Incidence light

Model	Factory setting
Unshown	Reflection light measurement E External incidence measurement

\*Important point

- ① May need angle calibration if you purchase in product-specific sale. What is angle calibration?⇒P.5-074
- ② We will calibrate at the factory shipping if you purchase it with autocollimator. What is set?⇒P.5-024

High speed·  
High resolution

Technical guide

Selection guide

Compact

Standard

High resolution wide view

Screen

Ghost cancel

2 wavelength

Processing unit

Accessories

## Processing Unit : HIP-1200

Angle measurement

Wavefront sensor

6-axis sensor O-PIAS

Laser autocollimator

Laser autocollimator set

Accessories for alignment



### Features



### ■Multi spots measurement

Several spots are analysed by multi spot reflection at the same time.

### ■3 kinds of applications

Runout measurement:Analys the tolerance range of tilt and runout.


3 beam measurement:For optical pick up grating.

Off set tilt:This function is useful when measuring the initial angle followed by relative angle in reference to the initial position.

### ■Automation

I/O for control and serial inter face are equiped as standard.

### ■Application




### ■Tilt measurement of image sensor for cameras.

Reflection of cover glass surface and element surface are captured at the same time.

Distinguish the sizes are sorted from the smallest shows tilt of element surface.

It is possible to judgment of acceptance for standard values.(Element surface can not take a reflection or it might be taking other reflection depending on the imager.)

Optical pick up




### ■Measure sub-beam ( $\pm$ primary light) of optical pick up

Detect a diffraction light by grating, and then measure the angle different to the direction of rotation.

It is possible to judgment of acceptance for standard values.

Precision Motors/Optical disc



### ■Measure runout and tilt of precision motors and optical disc.

Runout and tilt is measured in the axial trace of reflection.

It is possible to judgment of acceptance for standard values.

SPEC			
Model <sup>※1</sup>	HIP-1200/***-P	HIP-1200/***	HIP-1200/***-V
Basic specifications	Compatible autocollimator	H350	H400/H600
	Image output	NTSC signal (BNC)	
	Serial I/F	RS-232C (D-sub 9P male)	
	Parallel I/F	Insulated type I/O, in:3port, out:1port	
	Control	Membrane switch • Serial command control • I/O control	
	Power source (Included AC adaptor)	AC100~240V±10% 50/60Hz / DC12V 1A or less	
	Size, Weight	160 (W) × 50 (H) × 105 (D) mm, 600g	
Accessory	Communication speed	5Hz	
	Connecting cable	Nothing <sup>※2</sup>	HBNC-2
	AC adaptor	HDC12V-1000MA	HDC12V-1000MA

※1 See CP.5-040 for compatible type.


※2 Connect with H350 attachment cable.

## Function

Model	HIP-1200		HIP-1000*
	The latest model	End-of-life model	
Angle measurement	Single measurement	○	×
	Multi measurement	○	○
	(Level ordering)	Sizing/angle size from the center of measurement	Labeling from upper left
Application function	Off-set tilt measurement	○	×
	3 beam measurement	○	×
	Runout measurement	○	○
Angle analysis	Center of area	○	○
	Luminance centroid	○	○
	Luminance peak	○	×
User-friendly function	Display the luminance value	○	×
	Zoom function	○	×
	Image rotation (XY axis switching)	○	×
	Reversed a mark	○	×
	Averaging	○	×
Admission decision	Circle	○	○
	Square (Rad/Tan)	○	○
Admission decision Judgment condition	Offset judgment	○	×
	Multi measurement Selected points	○	×
	Multi measurement All points	○	—
Serial communication	Command control	○	—

※HIP-500/HIP-1000 has stopped being produced. Please contact us for maintenance.

## The display (Multi relative measurement MULTIR.)




## Processing Unit : HIP-1200

### A wide range of centroid analysis.

There are 3 kinds of centroid analysis at when measure the angle. (Center of area, Luminance Centroid, Luminance Peak)

Center of area measurement, ideal for specular reflection such as mirrors or glasses. Luminance centroid measurement is good for measuring skew angle such as lens edge surface. Also HIP-1200 has measured luminance peak.




### Single measurement

Measuring only one spot reflection.

Display the angle it calculated all of luminance information on the screen when the several refractors is measured with a single measurement on the luminance centroid mode. (A single measurement that is center of area of specified label.)


Single measurement SINGLE

One spot



Single measurement SINGLE

Several spots



(When selecte luminance centroid mode)  
Angle display a centroid that included all of luminance  
on the screen


### Multi measurement

Selectable from 2 kinds of measurements depending on your opinion.

MULTIA: Several spots are detected at the same time.


MULTIR: Angle difference of several spots are detected.

Multi measurement MULTI A



Measure an angle of several spots at the same time.

Multi relative measurement MULTI R



Measure an angle difference between spots.


## ■Labeling function

The arrangement of the label number is set.

Spots are numbered starting with the spot having the largest area (in the descending order of area).


Spots are numbered starting with the spot having the largest angle from the center of measurement (in the ascending order of angle).

### ■Sizing



Detect a size of spot, and make a label to be ordered from largest.

### ■Order from smallest




Spot is numbered 1 from smallest(near center of measurement). Label should be changed when spot is switched.


Centroid analysis Mode	Single measurement		Multi measurement		
	Setting code	Labeling	Centroid analysis	Labeling	Centroid analysis
Center of area	BIN	Label	Center of area for a specified label or Each center of area for all the labels	Label	Center of area for a specified label or Each center of area for all the labels
Luminance centroid	GRAY	Unlabel	Luminance centroid of all of the luminance in the screen	Label	Luminance centroid of one point of specified label or Each luminance centroid for all the labels.
Luminance peak	PEAK	Unlabel	Max. luminance position of all of the luminance on the screen.		It can not be used

## ■User-Friendly function

### ■Function 1 Rotation




Rotate the image without autocollimator head condition. ( $\pm 90$  deg.)




Freely attachable the image.  
Reduce your strange feeling easily.


### ■Function 2 Mirroring



Reverte the image without autocollimator head condition.




### ■Function 3 Offset judgment



Can be offset the judgment range to any position.

### ■Function 4 Zoom



Zoom magnification allows, and the display is zoomed automatically when the spot enters the zoom area.

### ■Function 5 External interface



I/O and command controllable. Control such as reset and data request are available from external equipment.

## Laser Autocollimator Set Type



### ■Features

#### ■6 kinds of setting types

Allows combination of multiple units for your need.

#### ■Included everything what you want.

Set model is included all what you need.

#### ■Accurate calibration

We are trying sensitive calibration for each our set model.

### ■Application



#### 【Full set】

See page ▶ P.5-051 for details.

Autocolli

HIP

Monitor

Monitor  
bracket

Tilt

Standard  
stand



#### 【High rigidity set】

See page ▶ P.5-053 for details.

Autocolli

HIP

Monitor

Monitor  
bracket

Tilt

High rigidity  
stand



#### 【monitor&HIP set】

See page ▶ P.5-055 for details.

Autocolli

HIP

Monitor

Monitor  
bracket

Tilt

Standard  
stand



#### 【HIP set】

See page ▶ P.5-057 for details.

Autocolli

HIP

Monitor

Monitor  
bracket

Tilt

Standard  
stand



#### 【Half set】

See page ▶ P.5-059 for details.

Autocolli

HIP

Monitor

Monitor  
bracket

Tilt

Standard  
stand



#### 【Monitor set】

See page ▶ P.5-061 for details.

Autocolli

HIP

Monitor

Monitor  
bracket

Tilt

Standard  
stand


※See autocollimator accessory if purchase separately. See ▶ P.5-063

Common accessories ▶ P.5-067


■ Vertical illumination position(e.g. Full-set or Half-set model)

※Please contact us about high rigidity model


■ H350R




■ H400-C




■ H400-CS




■ H350B



■ H450



■ H600



**DISCONTINUED**

■ Model

## H400-C025S/F12

1

2

1 Laser autocollimator model

Model	Product name
H350*-C****	Compact type
H400 -C****	General-purpose type
H400 -C****S	High resolution Two-visual-field
H450*-C****	Ghost cancel type
H600*-C****S	Dual-wavelength type

2 Set model

Model	Product name	Autocollimator	Processing unit	LCD monitor	Monitor bracket	Tilt stage	Stand	Accessary
/ F 12	Full-set 1200	○	HIP-1200	○	○	○	HA10	—
/ F 5 5	Full-set 550	○	HIP-550	○	○	○	HA10	—
/ K 1 2	High rigidity set 1200	○	HIP-1200	○	○	○	HA11N	—
/ K 5 5	High rigidity set 550	○	HIP-550	○	○	○	HA11N	—
/ M 1 2	Monitor and HIP set 1200	○	HIP-1200	○	—	—	—	—
/ M 5 5	Monitor and HIP set 550	○	HIP-550	○	—	—	—	—
/ H 1 2	HIP set 1200	○	HIP-1200	—	—	—	—	—
/ H 5 5	HIP set 550	○	HIP-550	—	—	—	—	—
/ D - H	Halfset (Degree display)	○	—	○	○	○	HA10	Degree displayed film
/ M - H	Halfset (Minute display)	○	—	○	○	○	HA10	Minute displayed film
/ D	Monitorset (Degree display)	○	—	○	—	—	—	Degree displayed film
/ M	Monitorset (Minute display)	○	—	○	—	—	—	Minute displayed film
Unshown	Laser autocollimator	○	—	—	—	—	—	—

# Optical Sensor

## Full Set

Autocollili    HIP    Monitor    Monitor bracket    Tilt    Standard stand



### ■ Features

#### ■ Select the standard model to help you decide

This model covers everything for measurement with autocollimator.

Satisfy multiple and diverse users' needs.

#### ■ Easy to use and accurate measurement


Tilt of work that is detected on the autocollimator by using processing unit, is displayed as numeric values on the XY two dimension coordinates.

#### ■ Two types processing unit

HIP-1200 is focused on the high performance, HIP-550 is placed especially for the price.

Selectable from both of our high grade processing units.

### ■ Switch • Conector • connection



※ See autocollimator accessory if purchase separately. See P.5-063  
Common accessories P.5-067



# Optical Sensor

## Monitor and HIP Set

Autocollil.    HIP    Monitor    Monitor bracket    Tilt    Standard stand



### ■ Features

#### ■ For embeded equipment

They seem to be a set; a laser autocollimator for measurement, a processing unit for analysis and then a LCD for display. You can also easily set up on production line or equipment integration with optional jig.

#### ■ Easily & Accurately


Display a numeric value on the XY, 2D coordinate workpiece tilt that is detected by laser autocollimator of processing unit. Can be judged angle result without dependence operator.

#### ■ 2 types of processing unit

Selectable from the HIP set 1200 focusing on the function and an emphasis cost performance the HIP set 550.

Both of 2 types have a communication interface as standard, correspond with the automation.

### ■ Switchi • Connector • Connection



※ See autocollimator accessory if purchase separately. See P.5-063  
Common accessories P.5-067



## Accessory for Autocollimator

Available accessories for each types such as beam split adaptor, 90degree folded adaptor or additional items.

Accessories index for each type		Screen type H100 series	Compact type H350 series	General-purpose type H400-C series	High rigidity two-visual-field type H400-CS series	Ghost cancel type H450 series	Dual-wavelength type H600 series	High speed high resolution type H900 series	Processing unit type HIP series
Pin-hole  ▶P.5-065			●	●	●	●	●		
Beam split 90 degree folded adaptor  ▶P.5-064			●	● <sup>*1</sup>	●	●	●		
ND filter  ▶P.5-065			Please contact us	Please contact us	Please contact us	●	●	Please contact us	
ND filter folder  ▶P.5-065			Please contact us	Please contact us	Please contact us	●	●	Please contact us	
Adaptor plate  ▶P.5-064				●					
Power source  ▶P.5-066			●	●	●	●	●	●	●
Image cable  ▶P.5-066				●	●	●	●	●	●
Connector, cable, mouse for processing unit  ▶P.5-066									●
Accessories for high rigidity H900  ▶P.5-066								●	

\*1 Extra plate for attachment of HT(for H350) is required. ▶P.5-064

Common accessories ▶ P.5-067	Tilt stage	Stand	Rod bracket
LCD monitor	Monitor bracket	Mirror	Wedge substrate


## Beam Split•90 Degree Folded Adaptor




90 degree folded adaptor



3 beam split adaptor



Folded 90 degree optical axis of autocollimator



Multi spot measurement is available with splitting an optical axis 2 or 3.


Product name	Model	Weight	Remark
90 degree folded adaptor	<b>HT-90</b>	0.04kg	Folded 90 degree optical axis of autocollimator. <sup>*</sup> 1 <sup>*</sup> 3
2 beam split adaptor	<b>HT-20</b>	0.24kg	Split two the autocollimator beam <sup>*</sup> 2 <sup>*</sup> 3 <sup>*</sup> 4
3 beam split adaptor	<b>HT-30</b>	0.25kg	Split three the autocollimator beam <sup>*</sup> 2 <sup>*</sup> 3 <sup>*</sup> 4

<sup>\*</sup>1 The folded direction can be changed as required.

<sup>\*</sup>2 Minimum split 14mm, maximum split 55.5mm between 2 beams

<sup>\*</sup>3 Extra plate for attachment of HT(for H350:350-HTP) is required.

<sup>\*</sup>4 Beam power variation occurs due to installation direction. (Contact our support team)



## Adaptor Plate

The plate for changing installation direction is used to mount the tilt stage on side of H350.

The plate for installation direction (H350-AP) is supplied with the full-set and half-set.

HT installation plate is used to mount the split or folded adaptor.




Plate for changing position(for H350 setting)


Product name	Model	Weight
Plate for changing position(for H350)	<b>H350-AP</b>	0.10kg
Plate for HT mounting(for H350)	<b>H350-HTP</b>	0.10kg

### [Mounted surface and display direction]


How to use the plate for changing position(for H350AP)



#### Set on the front (without H350-AP)




When mirror is tilted left ⇒ B




When mirror is tilted forward ⇒ C

#### Set on the side (with H350-AP)



When mirror is tilted left ⇒ C



When mirror is tilted forward ⇒ D



<sup>\*</sup> Autocollimator series other than H350 doesn't need an adaptor plate.

<sup>\*</sup> Note that H350 is different of installation and display direction other than autocollimator.

# Optical Sensor

## Pinhole

The pinhole for changing the outgoing beam of autocollimator. (As a standard). We provide  $\phi 0.5$  type<sup>※3</sup> of micro region measurement.



Product name	Model	Weight	Attached as standard <sup>※1</sup>
Pinhole for H100 ( $\phi 1, 2, 3, 4$ )	H100-P	0.02kg	H100
Pinhole for H350B/R-C100/C175 ( $\phi 0.5$ )	H350-C100-P05	0.01kg	
Pinhole for H350B/R-C050 ( $\phi 1, 2$ )	H350-C050-P12		H350 * -C050
Pinhole for H350B/R-C100/C175 ( $\phi 1, 2$ )	H350-C100-P12		H350 * -C100
Pinhole for H350B/R-C050 ( $\phi 3, 4$ )	H350-C050-P34		
Pinhole for H350B/R-C100/C175 ( $\phi 3, 4$ )	H350-C100-P34		
Pinhole for H400-C/H400-CS ( $\phi 1, 2, 3, 4$ )	H400-P	0.02kg	H400-C/H400-CS
Pinhole for H400-C (old type) ( $\phi 1, 2, 3, 4$ ) <sup>※2</sup>	H400-C-P		
Pinhole for H450 ( $\phi 1, 2, 5, 7$ )	H450-P		H450
Pinhole for H600 ( $\phi 1, 2, 3, 4$ )	H600-P		H600
Pinhole for H400-C/H400-CS ( $\phi 0.5$ )	H400-P05		
Pinhole for H450 ( $\phi 0.5$ )	H450-P05		
Pinhole for H600 ( $\phi 0.5$ )	H600-P05		

※1 Equipped it as standard at the factory shipping.

※2 H400-C upper connecting plate White(old type)/Black(new type)

※3  $\phi 0.5$  type : Not available H350B/R-C050.

## ND Filter



The ND filter for controlling light volume when measure external incident light.

ND filter (For H400-CS, H600)	
Size : $\phi 30$	
S73-30-1	Transmittance : 0.1% or less
S73-30-1P	Transmittance : 1%
S73-30-3	Transmittance : 3%
S73-30-13	Transmittance : 13%
S73-30-25	Transmittance : 25%
S73-30-50	Transmittance : 50%


※Not available reflection measurement in using internal light source

## ND Filter Holder



For H400-CS and H600-CS.

ND filter holder (For H400-CS/H600)	
H400-NDH	0.04kg
※Available $\phi 30$ ND	



## Electrical Power



HDC6V-1000MA



HDC12V-500MA



HDC12V-1000MA



HDC-CABLE1

Product name	Model	Weight	Remark
AC adapto r DC6V 1000mA	<b>HDC6V-1000MA</b>	0.08kg	For H350 (except H350R-C175)
AC adapto r DC12V 500mA	<b>HDC12V-500MA</b>	0.07kg	For H100, H400, H450, H600, H350R-C175
AC adapto r DC12V 1A	<b>HDC12V-1000MA</b>	0.10kg	For HIP-550, HIP-1200, ProCo, O-PIAS, VCM-562
DC power cable	<b>HDC-CABLE1</b>	0.01kg	Available power supply to LCD from HIP-550/HIP-1200

## Electrical Power



HBNC-1



HBNC-2



BNC-3C2V-2M

Product name	Model	Weight	Remark
Video cable for H350*-C050/C100	<b>HBNC-1</b>	0.06kg	BNC-POutput
Video cable for H400/H450/H600	<b>HBNC-2</b>	0.08kg	BNC-POutput
Video cable for H350R-C175	<b>HBNC-3</b>	0.06kg	BNC-POutput
BNC cable (2m)	<b>BNC-3C2V-2M</b>	0.15kg	BNC-P/BNC-P
Video cable for H350 (HIP500)	<b>H350-HIP</b>	0.10kg	Connector for HIP-500


## Connector Cable for Processing Unit



HIP-R03-PB5M



HIP-MOUSE



D100-R9-2

Product name	Model	Weight	Remark
I/O Connector	<b>HIP-R03-PB5M</b>	0.02kg	Compatible HIP-1200
Mouse for HIP-1000/2000/2400	<b>HIP-MOUSE</b>	0.15kg	D-Sub with a conversion connector
RS-232C Cross cable (9 Pin)	<b>D100-R9-2</b>	0.15kg	

## Accessories for High Rigidity H900(HRAD)



H900-CAP



H900-CABLE1



H900-CRING


Product name	Model	Weight	Remark
Cap for adjustment of optical axis	<b>H900-CAP</b>	0.01kg	Common H900
Sensor head cable (2m)	<b>H900-CABLE1</b>	0.22kg	Common H900
Ring for C mount	<b>H900-CRING</b>	0.01kg	For H900-P017

# Optical Sensor

## Configuration of The Common Accessory

This is the common accessory for laser autocollimator, ProCo and O-PIAS.

We have also laser autocollimator set model included the common accessory.⇒ Set format P.5-049




## Tilt Stage

RoHS

This is good for controlling of light axis tilt.  
This is useful for 2-axis tilt adjustment and lock mechanism.




SPEC			
	2-axis small-sized tilt stage	2-axis middle-sized tilt stage	2-axis big-sized tilt stage
Model	HB10	HB11	HB12
Load capacity	1kgf	2kgf	5kgf
Tilt angle	2 degrees	2 degrees	5 degrees
Material - finished	Aluminum-Almite finish, Stainless		
Weight	0.2kg	0.4kg	0.5kg



Compatible optical sensor	H100		○
	H350		○
	H400-C		○
	H400-CS	×	○
	H450	×	○
	H600	×	○
	H900 (HRAD)	×	○
Compatible stand	O-PIAS	×	○
	HA10	○	×
	HA10L	○	○
	HA11N	○	×
	HA12S	×	○
	HA12	×	○
	HA12L	×	○
Compatible rod bracket	HB10-RB20	○	×
	HB10-RB38	○	×
	HB12-RB20	×	○
	HB12-RB38	×	○

## Stand

Can be installed any optical sensors.  
It consists of stand base, pole and rod bracket.



SPEC						
	Stand (Aluminum base)	Stand (Aluminum long-base)	High rigidity stand (S45C-made base)	High rigidity stand (S45C-made long base)	High rigidity stand (S45C-made base)	Large High rigidity stand (S45C-made long base)
Model	HA10	HA10L	HA11N	HA12S	HA12	HA12L
Weight	2.4kg	3.0kg	9.2kg	9.2kg	11.4kg	19.1kg
Compatible tilt stage	HB10	○	○	○	×	
	HB11	○	○	○	×	
	HB12	×	○	×	○	

### Configuration

Paul	Model	A42-400		HA11-P38L355	HA12-P38L600
	Diameter	$\varphi 20$ (g6)		$\varphi 38$ (g6)	$\varphi 38$ (g6)
	Length	L=400mm		L=355mm	L=600mm
	Material	Stainless		Stainless	Stainless
Stand base	Model	HA10-BASE	HA10L-BASE	HA11N-BASE	
	Size	W200×D250	W200×D350	W200×D250	
	Material - finished	Aluminum-Almite finish		S45C electroless nickel plating	
Rod bracket	Model	HB10-RB20	HB10-RB20 (With plate)	HB10-RB38	HB12-RB38

## Rod Bracket

The rod bracket for installation tilt stage to the pole.

There are two kinds of radius  $\varphi 20$  and  $\varphi 38$ . Also available 2 kinds of tilt stages for HB10 and HB12.(A rod bracket is included in a stand)



SPEC				
	$\varphi 20$ Rod bracket	$\varphi 38$ Rod bracket	$\varphi 20$ Rod bracket	$\varphi 38$ Rod bracket
Model	HB10-RB20	HB10-RB38	HB12-RB20	HB12-RB38
Material - finished	Aluminum-Almite finish, Resin			
Weight	0.1kg	0.2kg	0.2kg	0.3kg
Compatible pole diameter	$\varphi 20$ (g6)	$\varphi 38$ (g6)	$\varphi 20$ (g6)	$\varphi 38$ (g6)
Compatible tilt stage	HB10	○	×	
	HB11	○	×	
	HB12	×	○	

# Optical Sensor

## LCD Monitor

RoHS

Angle  
measurement

Wavefront  
sensor

6-axis sensor  
O-PIAS

Laser  
autocollimator

Laser  
autocollimator  
set

Accessories  
for alignment

High speed·  
High resolution

Technical  
guide

Selection  
guide

Compact

Standard

High resolution  
wide view

Screen

Ghost  
cancel

2 wavelength

Processing  
unit


Accessories

5.6 inch LCD. Note that it is different model according to connect optical sensor directly and via processing unit.




(Include dedicated stand in LCD)


### Connect to H350 directly (visual observation-based)



### Connect to H400/H450/H600 directly (visual observation-based)



### Connect to HIP-550/HIP1200



### SPEC

	VCM-562P	VCM-562W	VCM-562HIP
Model			
Screen size		5.6 inch	
Pixel number		320 (W) × 3 (RGB) × 234 (H)	
power consumption		420mA(DC12V)、 5.0W	
Input system		AUTO (NTSC/PAL)	
External dimension		163 (W) × 112 (H) × 30 (D)	
Image cable length		1.5m	
Weight (except AC adaptor)		320g	

### Configuration

LCD monitor Mainbody	VCM-562-M		
BNC - RCA Conversion connector	BNC-P	BNC-W	
Power source	HDC12V-1000MA		HDC-CABLE1
Accessarry	Stand (self-supporting)		


Connect a monitor	H350	○		
	H400-C		○	
	H400-CS		○	
	H450		○	
	H600		○	
	HIP-550			○
	HIP-1200			○
	O-PIAS		○	

※The product and product specifications are subject to change without notice for modification.


## Monitor Bracket

Can be fixed the monitor to monitor bracket on  $\phi 20$  and  $\phi 38$  pole.

(Include monitor bracket in laser autocollimator full-set, half-set and high rigidity-set)



With stand



With monitor bracket

SPEC		
	$\phi 20$ Monitor bracket	$\phi 38$ Monitor bracket
Model	<b>HA10-MB</b>	<b>HA11-MB</b>
Material - finished	Aluminum-Almite finish、Resin	
Weight	0.3kg	0.4kg
Compatible pole diameter	$\phi 20$ (g6)	$\phi 38$ (g6)
HA10	○	×
HA10L	○	×
HA11N	×	○
HA12S	×	○
HA12	×	○
HA12L	×	○



## Electrical Power•Conecotor• Cable



HDC12V-1000MA



HDC-CABLE1

SPEC		
	AC adaptor DC12V1A	DC power source cable
Model	<b>HDC12V-1000MA</b>	<b>HDC-CABLE1</b>
Input	AC100~240V 50/60Hz	DC12V
Output	DC12V 1A	DC12V
Remark	For LCD(VCM-562)	Available supplying power to LCD from processing unit
Weight	0.10kg	0.01kg



BNC-3C2V-2M




BNC-W



BNC-P

SPEC			
	BNC Cable (2m)	BNC-RCA Conversion connector	
Model	<b>BNC-3C2V-2M</b>	<b>BNC-W</b>	<b>BNC-P</b>
Specification	BNC-P/BNC-P	RCA-P/BNC-P	RCA-P/BNC-J
Weight	0.15kg	0.05kg	0.05kg



## Mirror

For setting an optical sensor in vertically on the reference face.



SPEC	
	Parallel mirror
<b>Model</b>	<b>HS-0</b>
Material	BK7
External dimension	$\phi 30+0/-0.1\text{mm}$
Thickness	$10\pm0.1\text{mm}$
Parallelization	Less than 5 seconds
Surface accuracy	Duplex reflective wavefront $\lambda/10$
Polishing effective diameter	Center over $\phi 27\text{mm}$ circle
film forming specification	Dielectric multilayer film
Incident angle	$\theta=0^\circ$
Refrectance	$R650\pm20\text{nm}\geq99\%$
Film effective diameter	Center over $\phi 28\text{mm}$ circle

## Wedge Substrate/Wedge Mirror



Wedge substrate and Wedge mirror with micro angle.




Wedge millor

SPEC		
	Wedge substrate	Wedge mirror
Material	BK7	BK7
External dimension	$40\pm0.1\times40\pm0.1\text{mm}$	$40\pm0.1\times40\pm0.1\text{mm}$
Thickness	10mm	10mm
Wedge tolerance	$\pm10$ seconds	$\pm10$ seconds
Surface accuracy	Duplex reflective wavefront $\lambda/4$	Duplex reflective wavefront $\lambda/4$
Surface accuracy measurement wave length	$\lambda=632.8\text{nm}$	$\lambda=632.8\text{nm}$
Polishing effective diameter	Center over $\phi 36\text{mm}$ circle	Center over $\phi 36\text{mm}$ circle
Film forming specification	N/A	AL+MgF2

	Wedge substrate	Wedge mirror
Angle	Model	Model
1 degree	HS-100	HS-100AL
0.5 degrees	HS-050	HS-050AL
0.25 degrees	HS-025	HS-025AL
0.2 degrees	HS-020	HS-020AL
0.1 degrees	HS-010	HS-010AL




## Laser Autocollimator Technical Guide

### ■What is a Laser autocollimator?

#### Laser autocollimator


- Measure instantaneously the angle in a single location.
- $\theta_{xy}$  2 demension coordinates display
- Long width
- High resolution 8 sec.  
(Measurement range  $\pm 0.5$  degree)
- Changeable measurement beam radius.



Measuring object is only regular reflection of mirror and glass surface.  
See our catalog P5-026 for scatter reflections such as resin or metal.


#### Length measurement system

Calculates the differences in height between any 2 points with a length measurement sensor.



### ■Angle measurement principle

When the parallel light incents to the collecting lens, it will be converted to the position information on the focal plane depending on the incidence angle. Irradiation of parallel laser beam to the measurement object using this principle, the tilt value of the object will be measured by collecting reflection light on the CCD.



Left drawing: If the object is perpendicular to the irradiation beam, the reflection light returns to vertical, and light will be collected in the center.


Right drawing: In case of the object is tilted, the reflection light will tilt to the left, and also collecting area will move to the left.

Can be known the tilt by measuring the different between the reflection light and collecting area.

### ■Unchanged the measurement result, even the measurement distance is changed.

The measurement result that is gotten by using an autocollimator as an angle is measured directly even if the distance between the object is changed. Available for pitching and yawing measurement of the up and down drive unit.

Imaging in the same point




### ■Specialized usage: External light source incidence measurement

#### Normal (The reflection measurement)

Normally an autocollimator is used in the reflection measurement. (Left drawing) When the object tilt is  $\theta$  degree, the reflection light returns  $2\theta$  degree angle. Autocollimator shows tilt angle  $\theta$  degree of the object.

(e.g.) When the object tilts 0.1-degree, 0.2-degree tilted reflection light returns to an autocollimator, but measurement value is displayed 0.1-degree.



#### Measure the external light source

Can be measured laser beam tilt that is emitted from the object light source by turning off the semiconductor laser (LD) in the autocollimator.


## ■Display the angle

Display the spot position that is imaged on a screen and CCD on the 2D angular coordinates.

CCD built-in type (H350/H400/H450/H600) displays X · Y · D in digital by using the optional processing unit (monitor).

The D value shall be calculated as follows.

Selectable angular unit from degree, minute or second.




## ■Viewing range

Measurement range (viewing) of laser autocollimator is decided by the designed value of optics.

We have a variety of range  $\pm 0.16^\circ \sim \pm 2^\circ$  as belows

"H\*\*\*-C16" in the model number shows  $\pm 0.16^\circ$



## ■Pixel resolution

Pixel resolution is decided by the number of CCD pixels in CCD built-in type (H350/H400/H450/H600).

Visual field range  $\pm 0.5$  degree (e.g.H400-C50)

$\pm 0.5^\circ = 1^\circ$  (Range)


$1 \div 480 = 0.002083^\circ$

$= 0.125$  minute

$= 7.5$  seconds

$\doteq 8$  seconds

Below decimal point rounded




Available higher resolution measurement of the sub-pixel processing in the centroid analysis mode of processing unit.

# Laser Autocollimator Technical Guide [Notice]

## ■Object measurement distance

Using it move away from measurable area, measurement range cannot be covered whole area due to reflection through the outside of aperture.

Please use the product as specified within the data.



## ■About semiconductor laser life

LD that is internal semiconductor is consumable. Life guarantee is a year (8760 hours) under recommended environment.

# Optical Sensor

## Laser Autocollimator Technical Guide [Measurement Technique]

Angle measurement

Wavefront sensor

6-axis sensor O-PIAS


Laser autocollimator

Laser autocollimator set

Accessories for alignment


### ■ Measurement of right-angle prism

Absorb the vertical plane of the right-angle prism, and the measurement of the right angle measures the horizontal surface with laser autocollimator.




No limit setting direction of the main autocollimator.

Can be measured the vertical plane of the right-angle prism in the horizontal setting.




### ■ Measurement of scattering reflector (resin or metal)

Laser autocollimator cannot be measured directly due to a collectible the reflection on the scattering and curvature surface.




Can be measured the objective tilt with a parallel mirror on the measurement surface.



### ■ Simultaneous measurement of dihedral parallelism


It is possible to compare the angle even if the object height or position is different.



《Right》Measure parallelism between the dihedral on coaxial by transmitting a laser beam.

《Left》Measure parallelism of uneven dihedral by using 2 beam branch adaptor (optional).

Both reflection spot might be overlapped a point if it is parallelism.



### ■ Lens measurement

Lens is measured by focussing laser beam to the edge.


Problem:

Narrow the edge.

Not easy to measure due to the ghost refraction.


<Problem solving>

Optical system H450 (O-035) will be irradiated whole lens (max.  $\phi 8$ ) beam, then measures only flat surface that returns from the edge.



Problem :

The edge is rough and not mirror surface.




<Problem solving>

Changeable lens aperture by using the jig for receiving the edge of the lens.

## ■Centroid analysis tip

The luminance centroid mode is good for measuring the tilt of lower flatness such as the lens edge.



Compared with the area centroid mode, measurement linearity and reputability will be improved dramatically.




Use when the luminance centroid mode is within the appropriate range constantly checking the luminance bar for enhanced measurement accuracy.

## ■Runout measurement of precision motor and optical disc drive

Selectable from HIP-1200 (cost-sensitive) (C) 5-041 and HARD H900 (accuracy-oriented) (C) 5-017 for runout analyzing of precision motor and optical disc drive.



## ■Glossary: Binarization threshold




The standard for changing the light detected by the autocollimator optical system from grayscale image to a binary image.

Binarization processing is the image processing method. When the pixel is 1 (white) that over standard the luminance value, the other is 0 (black). This standard (Minimum blights for being white) is called Binarization.

※ Note: Set up with area centroid mode.

## ■Glossary: 256 steps



The luminance that is detected by the autocollimator optical system will be displayed in 8 bits grayscale image.

The grayscale image centroid is calculated by luminance centroid analysis.


Compared with the area centroid analysis with binarization processing, can be analyzed accurately.

※ Note: Threshold must be set in the luminance centroid mode due to avoid basing noise.

## ■Optical disc drive assembling

This makes it easier to control the optical disc drive in HIP-1200.

(C) 5-041 Drive the reflection of main or sub axis at the center of tilt of the motor runout measurement. Switchable easier runout and angle measurement by using I/O control. Reflect an offset skew value is controlled by a serial command.



High speed·  
High resolution

Technical  
guide

Selection  
guide

Compact

Standard

High resolution  
wide view

Screen

Ghost  
cancel

2 wavelength

Processing  
unit

Accessories