

RIMA™ RAMAN IMAGING SYSTEM

MEGAPIXEL IMAGES IN MINUTES!



The perfect Raman imager for the analysis of nanomaterials from graphene to carbon nanotubes, RIMA is a state-of-the-art ultrafast hyperspectral imaging system available at various excitation wavelengths (532 nm, 660 nm, 785 nm). RIMA is also a tool of choice for non-invasive monitoring and analysis of biological tissue.



RIMA NANO - 532 nm, 660 nm

TECHNICAL SPECIFICATIONS			
	RIMA 532	RIMA 660	RIMA 785
Spectral Range*	190 to 4000 cm^{-1}	100 to 4000 cm^{-1}	130 to 3200 cm^{-1}
Spectral Resolution	< 7 cm^{-1}	< 6 cm^{-1}	< 5 cm^{-1}
Microscope	Upright	Upright	Inverted
Objectives	20X, 50X, 100X	20X, 50X, 100X	20X, 60X, 100X
Excitation Wavelengths*	532 nm	660 nm	785 nm
Spatial Resolution	Sub-micron		
Maximum Scanning Speed	250 $\mu\text{m}^2/\text{min}$ at full spectral range		
Wavelength Absolute Accuracy	1 cm^{-1}		
Camera*	Back-illuminated CCD or sCMOS camera 1024x1024 px		
Video Mode	Megapixel camera for sample visualization		
Preprocessing	Spatial filtering, statistical tools, spectrum extraction, data normalization, spectral calibration		
Hyperspectral Data Format	FITS, HDF5		
Single Image Data Format	JPG, PNG, TIFF, CSV, PDF, SGV		
Software	Computer with PHySpec™ control and analysis software included		

UPGRADES*	RIMA 532	RIMA 660	RIMA 785
	Low-Noise Back-Illuminated Camera, EMCCD	Low-Noise Back-Illuminated Camera, EMCCD	Deep-depletion camera, EMCCD
	Additional excitation wavelengths available	Additional excitation wavelengths available	Broadband COL Camera, Motorized stage with piezo positioning on z-axis
	Spectral Range Extension: Anti Stokes	Broadband COL Camera: Color 3MP Camera	Spectral Range Extension: Anti Stokes
	Broadband COL Camera: Color 3MP Camera	FIGURE 1	Additional excitation wavelengths available

Hyperspectral Raman imaging using Bragg tunable filters of graphene and other low dimensional materials

Etienne Gaufres, Stéphane Marquet, Vincent Aymong, Nathalie Y-Wa Tang, Alexandre Favron, Felix Thouin, Charlotte Allard, David Rioux, Nicolas Cottenye, Marc Verhaegen and Richard Martel.

Journal of
**RAMAN
SPECTROSCOPY**




Figure A. (a) $130 \mu\text{m} \times 130 \mu\text{m}$ Raman mappings of the G peak intensity at $\lambda = 532 \text{ nm}$ of graphene bilayer islands on a graphene monolayer. (b,c) Spectra of monolayer (blue) graphene and of nonresonant (green) and resonant (red) bilayer graphene islands from selected points in (a). The peak indicated by * is an instrument artifact. (d) Raman image ($70 \times 47 \mu\text{m}^2$) of the G peak intensity of an artificial bilayer of graphene composed of two monolayers stacked on top of each other.




Figure B. (a) Raman spectrum at $\lambda_{\text{exc}} = 532 \text{ nm}$ of few layers MoS_2 extracted from a RIMA hyperspectral cube of the sample and corresponding to the area pointed by a cross in (b). (b) Color coded cartography ($130 \mu\text{m} \times 130 \mu\text{m}$) of the layer composition of exfoliated MoS_2 deposited on $100 \text{ nm SiO}_2/\text{Si}$ substrate. The color code is obtained from the difference in peak positions between the A_{1g} and E^{1}_{2g} modes.

Figure C. (a) $260 \times 260 \mu\text{m}^2$ Raman mapping of 6T molecules encapsulated in carbon nanotubes (6T@SWCNTs). (b) A representative Raman spectrum of the sample showing the characteristic peaks of 6T around 1460 cm^{-1} and the G band of CNTs around 1590 cm^{-1} .


Figure C. (a) $260 \times 260 \mu\text{m}^2$ Raman mapping of 6T molecules encapsulated in carbon nanotubes (6T@SWCNTs). (b) A representative Raman spectrum of the sample showing the characteristic peaks of 6T around 1460 cm^{-1} and the G band of CNTs around 1590 cm^{-1} .

Figure C. (a) $260 \times 260 \mu\text{m}^2$ Raman mapping of 6T molecules encapsulated in carbon nanotubes (6T@SWCNTs). The image is a superposition of the maximum intensity of CNTs at 1590 cm^{-1} (green scale) and 6T at 1450 cm^{-1} (red scale) obtained after background subtraction. Empty CNTs in green can be distinguished from filled CNTs with 6T molecules in yellow or red, depending on the intensity. (b) A representative Raman spectrum of the sample showing the characteristic peaks of 6T around 1460 cm^{-1} and the G band of CNTs around 1590 cm^{-1} . Adapted from [37].



Electrostatic Deposition of Large-Surface Graphene

Charles Trudeau, Laura-Isabelle Dion-Bertrand, Sankha Mukherjee, Richard Martel and Sylvain G. Cloutier




(a) White-light hyperspectral image with high field-of-view showing the edge of the deposition (dashed line). (b) Hyperspectral image of the full graphene deposition mapping the position of the highest intensity around the G peak ($1500\text{--}1600 \text{ cm}^{-1}$). The white box represents $130 \mu\text{m} \times 130 \mu\text{m}$. Acquired using RIMATM NANO - Photon Etc

Giant Raman scattering from J-aggregated dyes inside carbon nanotubes for multispectral imaging

nature
photonics

E. Gaufres, N. Y.-Wa Tang, F. Lapointe, J. Cabana, M.-A. Nadon, N. Cottenye, F. Raymond, T. Szkopek and R. Martel



Raman multiplexing, protein recognition and tagged bacteria with dyes@SWNTs nanoprobes. (a) Raman hyperspectral image at $1/4532 \text{ nm}$ of isolated bundles of 6T@SWNTs (red) and bcar@SWNTs (green) co-deposited at low coverage onto a Si/SiO_2 substrate. (b) As in a, but using a mixture of 6T@SWNTs, bcar@SWNTs and Ph@SWNT (blue) nanoprobes on Si/SiO_2 . (c) Top image: optical image of *Candida albicans* tagged with bcar@PEG-SWNT. Bottom image: corresponding Raman image taken at 532 nm of the bcar@f-SWNT mode centred at $1,520 \text{ cm}^{-1}$. (d) Raman image of the bcar@PEG-biot-SWNT probe taken at 532 nm using the peak centred at $1,520 \text{ cm}^{-1}$. The bcar@PEG-biot-SWNT probes selectively attached to immobilized streptavidin by microcontact printing in circular dot shapes (diameter, 10 nm). Inset: results using the reverse pattern with surface streptavidin located surrounding the dots.