# TUNABLE LIGHT SOURCES TLS-55-X300 TLS-72-X300 TLS-55-X300-SS TLS-55-Q250 TLS-72-Q250 TLS-55-X300-SS



# **Features**

- Monochromatic light from 300 nm up to 1800 nm
- Collimated, condensed, or coupled output light options. Light can be coupled to any optical system (including fibers and fiber bundles)
- Adjustable optical resolution from 20 nm down to 0.2 nm
- Flexibility of optional features: full automation, various sources and monochromators and more

1450 Global Drive, London, Ontario Canada, N6N 1R3 Phone: 519 644 0135 / Fax: 519 644 0136 Email: <u>sales@sciencetech-inc.com</u> www.sciencetech-inc.com

# **Applications**

- Absorption/transmission/reflection measurement systems
- As an excitation light source in fluorescence measurement systems
- Eye protection products
   measurement system
- PEC photochemistry measurement systems
- Solar cell quantum efficiency measurement systems



## I. Overview


Sciencetech offers a series of computer-controlled adjustable sources of monochromatic light, our Tunable Light Source (TLS) line.

Each TLS system is pre-aligned at the Sciencetech test lab and packaged with the system's test data and detailed manual. Sciencetech's tunable light sources are motorized instruments controlled by Sciencetech's TLS software. LabView drivers and ActiveX and DotNet modules are also available making the TLS a flexible system for integration with your experiment or Each base model comes with standard features:

- \* I" collimated light output
- \* Monochromator with two gratings
- \* Motorized filter wheel with 4 long pass filters installed
- \* Adjustable slits for optical bandpass control.
- \* Adjustable computer controlled power supply (Xenon lamps only)

\* SciSpec software for controlling the monochromator and filter wheel.

\* Reinforced optical base plate



instrument.

Sciencetech TLS systems are built with modular components making them flexible and adaptable well beyond the standard models offered.

| Product<br>Line | Monochromator | Lamp Type | Model       |
|-----------------|---------------|-----------|-------------|
| TLS             | -72 (1/8 m)   | -300×     | TLS-72-300X |
|                 |               | -250Q     | TLS-72-250Q |
|                 | -55 (1/4m)    | -300×     | TLS-55-300X |
|                 |               | -250Q     | LH-55-250Q  |

Choose –X light source for Xenon—best for UV/VIS Choose—Q light source for QTH—best for VIS/NIR There are several available output options: different focusing and collimating beams, different coupling with fibers/bundles, or coupling with different optical devices.

Available optional upgrades include: Motorized or manual Iris, motorized slits (manual slit is standard), Fixed slits, optical light feedback for additional source stabilization\*\*

#### **Customized Solutions**

All components of the system, including the light source and monochromator, can be changed. Available light sources include continuous

- Xe, Xe-Hg and Hg lamps from 75W-1.6 kW
- QTH lamps 50 W-500 W
- Dual deuterium/QTH source
- Mid-IR sources up to 1700Kelvin
- Available spectral range 180 nm 10 um \*
- Double subtractive and double additive monochromators
  - Monochromators up to 2.0 meter focal length



\* Light source, Optics and Grating dependent \*\* Only available with Xe and Xe-Hg Lamps

# 2. Specifications

Sciencetech's software, Sci-Spec, controls all components of the system. In the standard configuration, it controls the power supply of the light source, shutter, filter wheel, and monochromator. As an option, the user can add computer control on the input and output slits, and/or control of the iris. In the standard configuration, the output beam can be collimated or focused. Coupling with different devices is available as an option, such as coupling with a sample chamber, fiber, or fiber bundle.

|                              | TLS-72-X300  | TLS-72-Q250                                   | TLS-55-X300                        | TLS-55-Q250                 |
|------------------------------|--|---|------------------------------------|-----------------------------|
| Lamp Type                    | Xe 300 W   | QTH 250 W                                     | Xe 300 W                           | QTH 250 W                   |
| Monochromator Type           | 9072S (1/  | -   | 9055 (1/                           |                             |
|                              | , , , , , , , , , , , , , , , , , , ,  | ,   |                                    |                             |
| Functional Spectral range    | 300-1800 nm  |   |                                    |                             |
|                              | 0.5 nm @ 300-700 nm 0.2 nm @ 300-700 nm  |   | 300-700 nm                         |                             |
| Optical resolution           | 0.7 nm @ 700-1800 nm   |   | 0.4 nm @ 700-1800 nm               |                             |
| Filter wheel with filter set | Computer-controlled  |   | 6-position filter wheel            |                             |
| Shutter                      | Optional   | Optional computer-controlled shutter and expo |                                    | re control                  |
| Beam output                  | l" diameter collimated   |   |                                    |                             |
| Wavelength Repeatability     | 0.1 nm   |   | 0.03 nm                            |                             |
| Wavelength Accuracy          |  | 0.2   | nm                                 |                             |
| Intensity Control            | Source intensity adjustable  |   |                                    |                             |
| Intensity Control            | Optional manual iris   |   |                                    |                             |
|                              | 2 plane ruled gratings 30x30 mm  |   | 2 plane ruled gratings 50x50 mm    |                             |
| Gratings                     | 1200 gr/mm@300nm   |   | 1200 gr/mm@300nm                   |                             |
|                              | 600 gr/mm@1000nm   |   | 600 gr/mm@1000nm                   |                             |
| Bandpass                     | Two manual bilaterally-adjustable slits with vertical curtain attenuators                                    |   |                                    |                             |
| Optical Height               |  |   |                                    |                             |
| Power supply                 | Touchscreen, Con-<br>stant Current   | Manual, Constant<br>Current                   | Touchscreen, Con-<br>stant Current | Manual, Constant<br>Current |
| Software                     |  |   |                                    |                             |
|                              | USB—monochromator and filter wheel, RS232 Power Supply (xenon lamp only)                                     |   |                                    |                             |
|                              | IEC 61010-1:2010 Safety Requirements for electrical equipment for measurement,                               |   |                                    |                             |
| CE Certification             | control and laboratory use—Part 1  |   |                                    |                             |
|                              | IEC 61326-1:2012 Electrical equipment for measurement, control and laboratory<br>use—EMC requirements—Part 1 |   |                                    |                             |
|                              |  | use—EMC requ                                  | irements—Part I                    |                             |



### 3. Tunable Light Source—Solar Simulator

Sciencetech also offers a tunable light source solar simulator combination instrument. This combination system allows a selectable output at the target plane of spatially uniform sun light or monochromatic light.

The tunable light source—solar simulator combination system adds a homogenization unit to the exit port of the monochromator that includes and electronic shutter and 2 position filter tray.


The tunable light source functions the same as the other models in Sciencetech's tunable light source line but adds the extra functionality of a solar simulator with broadband spatially uniform white light.


Sciencetech currently offers one model of tunable light sourcesolar simulator:

#### TLS-55-X300-SS

The standard configuration of a TLS-SS includes a lamp housing with a elliptical reflector and xenon lamp, power supply with touchscreen, I/4 meter monochromator, automated filter wheel (with order sorting filters) and electronic shutter mounted to a metal breadboard to create a compact, fully assembled illuminator, controlled by original software. Each system is pre-aligned during production and packaged with the system's test data and detailed manual.


#### **Homogenization Unit**





#### **Solar Simulator Mode**

#### Tunable Light Source Mode






# 4. Specifications

Sciencetech's software, Sci-Spec, controls all components of the system. In the standard configuration, it controls the power supply of the light source, shutter, filter wheel, and monochromator. As an option, the user can add computer control on the input and output slits, and/or control of the iris. In the standard configuration, the output beam can be collimated or focused. Coupling with different devices is available as an option, such as coupling with a sample chamber, fiber, or fiber bundle.

| Model                        |  |  |
|------------------------------|--|--|
| Model                        | TLS-55-X300-SS   |  |
| Lamp Type                    |  |  |
| Monochromator Type           |  |  |
| Functional Spectral range    | 300-1800 nm  |  |
| Optical resolution           | Up to 0.2 nm @ 200-700 nm  |  |
|                              | Up to 0.4 nm @ 700-2500 nm   |  |
| Filter wheel with filter set | Computer-controlled 6-position filter wheel                                    |  |
| Shutter                      | Computer-controlled shutter and exposure control                               |  |
| Beam output                  | Spatially uniform output over 1" diameter                                      |  |
| Solar Simulator, Spectral    | Class A, AM1.5G  |  |
| Match                        | Class B, AM0   |  |
| Solar Simulator, Non-        | Class B  |  |
| Uniformity                   | Class B  |  |
| Solar Simulator, Temporal    | Class A  |  |
| Instability                  |  |  |
| Working Distance             | 80 to 140 mm   |  |
| Wavelength Repeatability     | 0.03 nm  |  |
| Wavelength Accuracy          | 0.2 nm   |  |
|                              | 2 plane ruled gratings 50x50 mm  |  |
| Gratings                     | 1200 gr/mm@300nm   |  |
|                              | 600 gr/mm@1000nm   |  |
| Slits                        | Two manual bilaterally-adjustable slits  |  |
| Optical Height               | 216mm (adjustable +/- 5mm)   |  |
| Power cumpl                  | Touchscreen,   |  |
| Power supply                 | Constant Current   |  |
| Software                     | Sci-Spec   |  |
| Computer Interface           | USB  |  |
|                              | IEC 61010-1:2010 Safety Requirements for electrical equipment for measurement, |  |
|                              | control and laboratory use—Part I  |  |
| CE Certification             | IEC 61326-1:2012 Electrical equipment for measurement, control and laboratory  |  |
|                              | use—EMC requirements—Part I  |  |



# 4. Specifications



# 5. Accessories

#### **Host Computer**

A host computer is required to operate Sciencetech's line of Tunable Light Sources. Software can be provided for installation on a customer's own host computer, or a preconfigured host computer system can be purchased at the time order, which comes with all software preinstalled and tested for easiest setup.

| SKU      | Product                                    |
|----------|--|
| 490-0128 | Preconfigured Host Desktop Computer System |
| 490-0127 | Preconfigured Host Tablet Computer System  |

### **Optional Accessories**

A number of accessories are available to upgrade the mechanical components of the system or add additional features.

| SKU                          | Product   |
|------------------------------|---|
| 810-0010                     | (TLS-CAL-PACKAGE) This service provides light source characterization<br>* power from 300-1800nm<br>* beam divergence<br>* light source stability ) |
| 120-9053                     | (SS-80-M) Motorized slit(s)   |
| Please inquire               | Manual or Motorized Iris  |
| 115-9027                     | (FS-02-N) Light Intensity Stabilizer  |
| 585-0002<br>With<br>585-0176 | Broadband Thermopile Detector, (BBT-03)<br>with<br>(UNO-1) Handheld Monitor for Thermopile Detectors  |


## **Output Accessories**

Different output accessories can be fitted to the monochromator output to customize the system for different applications—either as a stand-alone unit or to be used in conjunction with other components.

| SKU            | Product   |
|----------------|---|
| Please inquire | Focused Light Output (F/)   |
| Please inquire | Coupling to Fiber Cable/Bundle  |
| Please inquire | Optical Fiber (2m standard, custom lengths available); Optical bundle (custom length) |
| Please inquire | Custom requirements   |




# 6. Dimensions (mm)



## 9055 Monochromator-Based:

*Left*: The dimensions of the Sciencetech Tunable Light Sources are different for different models, but the 9055 monochromator-based model has the following dimensions in millimeters.



1450 Global Drive, London, Ontario Canada, N6N 1R3 Phone: 519 644 0135 Fax: 519 644 0136 Email: sales@sciencetech-inc.com www.sciencetech-inc.com

#### 9072 Monochromator-Based:

Left: The dimensions of the Sciencetech Tunable Light Sources are different for different models, but the 9072 monochromator-based model has the following dimensions in millimeters.

For TLS-55-X300-SS dimensions please contact a Sciencetech Application Scientist.

