Lower Back Laser

Nasal laser | Handheld Laser | 1x4 Laser | Neck laser | Neck laser 2.0
Leg Laser | Lower Back Laser | Wrist Laser

Lower Back Laser Wellness Device

The Waist Laser Therapy Device can be secured to the lower back, waist, or abdominal area using the included adjustable strap.

It uses 650nm + 790nm red light, delivered through laser diodes, to provide gentle massage and help relieve pain by directly delivering energy to areas of discomfort.

The treatment time can be freely set from 10 to 60 minutes, using continuous wave mode to promote blood circulation and alleviate discomfort.

Designed for continuous use, it's ideal for daily wellness and preventive care.



Key Features

Your go-to relief for lower back, pelvic, and abdominal discomfort

- Soothing warmth and gentle massage
- Eases lower back pain
- Relieves uterine and abdominal soreness
- Soothes pelvic and lumbar tension



Specification

Wavelength: 650 ±20nm / 790 ±20nm	Number of Laser Diodes: 98pcs
Product Size: L210/W130/H30 ±1mm	Laser Safety Classification:
Product Weight: 462 ±10g	Operating mode: continuous wave (CW)

Caution

- Avoid direct exposure of the laser output to the eyes to protect your vision.
- Do not power the device directly using a transformer, as this may cause electric shock or damage the main unit.
- Do not use fast charging adapters or operate the device while charging, as this may damage the unit.
- Keep the laser output panel clean during use to ensure effective energy delivery.
- Store the device out of reach of children and pets to prevent accidental ingestion or improper handling.
- If you are sensitive to light, have special medical conditions, or experience discomfort during use, stop using the product immediately and consult a physician.
- Do not use this device in combination with electronic medical equipment, as it may cause malfunction or safety hazards.
- Store in a cool, dry place, away from moisture and water. Handle with care to prevent damage from drops or impacts.
- If the product malfunctions or requires service, please contact customer support. Do not attempt to disassemble the device yourself.
- After charging is complete, or when not charging, disconnect the USB charging cable from the device.
- Keep the device away from microwaves, induction cookers, X-ray machines, CT scanners, and other strong electromagnetic equipment during use.
- Do not clean the product with solvents, thinners, gasoline, or similar substances, as this may cause malfunction or damage.
- Do not sterilize any accessories by boiling, chemicals, or microwave steam. Also avoid chlorine-based disinfectants or bleach.
- Handle the device with care during use. Do not drop, knock, crush, or apply strong impact to avoid damaging the unit.
- Do not use the device while it is charging.
- This product should be used at least once every 3 months. If not in regular use, it is recommended to charge the device once every 3 months using the standard charging method.

BACK

Need technical support or help?

CONTACT

Copyright © 2023 All Rights Reserved by LECC