

OlisWorks software for the HP 8453 Diode Array Spectrophotometer


Replace ChemStation with OlisWorks and enjoy all the benefits of a fast, easy to use data acquisition program in a Windows® 7/10 platform. Data export directly into Excel® or Olis GlobalWorks for 2D & 3D data processing.

Collection Modes Supported:

Single Wavelength Assay
Single Point Mode
Multiple Wavelength Assay
Single Scan
Repeated Scans
Continuous Scans
Rapid Scanning


Technical Specifications remain:

- Single beam
- 1 point per nm
- Deuterium & Tungsten lamps
- 190 – 1100 nm
- 1 nm resolution
- Faster rate of 10 scans per second over a span up to 100 nm



Benefits:

Gives the aged HP 8453 a new, simple to use interface
Utilizes a LAN connection for faster data transfer
No hardware changes are made to the 8453
Continued support by OLIS
Future software upgrades are optional, not mandatory



sales@olisweb.com or 706-353-6547