

VERE


Toll Free in USA
1-800-568-8373

[Home](#)[Optical Tables](#)[Instruction Manual](#)[Table Enclosures](#)[Ordering Information](#)[Price Quotation](#)


Over Table Shelf Assemblies

The Over Table Shelf Assembly is ideal for optimizing your otherwise wasted laboratory space. The 24" deep recessed shelf has a 2" grid of 1/4-20 holes on the underside for attachment of hanging shelves or a lighting system. Includes a power strip with 3-prong 110vac electrical outlets. The leveling feet are positioned under the optical table to allow un-obstructed movement around the table.


Part number: **OTSA24xx** (xx = your table length in inches)



Shown with optional hanging shelves



Top shelf as seen from underneath.
Observe grid pattern



Optional Accessories

Hanging Shelf

The hanging shelf is 12" high, 22" wide and 20" deep. Supports light loads up to 50 Lbs. Shelf is supported from underside by our support post (included).

Part number: **OTSA-HS-1000**




Optional shelf shown

Second Full Size Shelf

Positioned below top shelf. Contact our sales department for more details.

Part number: **OTSA-SS-xx** (xx = your OTSA length in inches)



Lighting Systems

Slim line under cabinet style light, 42"L x 5.5" W x 1" thick with (2)-T5 lamps

Part number: **OTSA-LSL**

Additional Power Strip

Part number: **OTSA-PS**



Optional Rear Support Version of the Over Table Shelf Assembly





3 Schreiber Industrial Park - New Kensington PA 15068 - 4549

Toll Free: 1-800-568-8373

Phone (724) - 335-5530 FAX (724) -335 -5569 Email Inquiries : [SALESdesk@ vere.com](mailto:SALESdesk@vere.com)