

orchestra[™] FROG

The *orchestra*[™] FROG (Frequency-Resolved Optical Gating) combines SHG (second harmonic generation) and TG (transient grating) FROG in one device with total spectral range from 200 nm to 2 μ m.

Wavelength range:	200nm to 2 μ m
TG range:	200nm - 1050nm
SHG range:	700nm - 2 μ m
Shortest pulses:	2 cycles from VIS to IR
High accuracy:	20 attosecond motor step
Long scan range:	200 picosecond
Polarization:	s or p
Small footprint:	16cm x 16cm

User oriented software for phase retrieval and pulse reconstruction; data compatibility with other software.



Address

1650 Boul. Lionel-Boulet
J3X 1P7 Varennes, Quebec
Canada

Copyright © 2021 *few-cycle*. All rights reserved.